

E.A.S.T. HP®

NEW!

MINOR SURGERY WITHOUT SCALES

A newly-conceived, practical and handy tool without wires, to be combined with the treatments performed with other devices present in the specialized studies of Dermatology, Plastic Surgery and Aesthetic Medicine.



Correction of palpebral ptosis of medium and slight...?

Soft Lifting of excess skin

Cutaneous hyperchromia

E.A.S.T. HP is widely used for phototype 4 and 5 (dark skin) without color problems

Striae periumbilical expanses

Warts, nevi, keloids, fibroids, xanthelomas, dyskeratosis

Correction of scars

E.A.S.T. (Electrical Arc Sublimation Therapy) is recommended in the outpatient treatment of excellence to eliminate skin imperfections in a secure and comfortable manner

EFFETS AND BENEFITS

CLASSE IIB
DIR. 93/42/CEE
e S.M.I. (Medical
Devices) CE 0068

THE ONLY DEVICE
IN THE WORLD WITH A DUAL
CYCLE SYSTEM WITH BATTERY
MICROPROCESSOR

FEATURES

- ▶ Dermal sublimation micro radio scalpel with voltaic arc
- ▶ It's practical and ergonomic
- ▶ Rechargeable lithium ion batteries with a overcurrent and overcharge circuit of high Italian quality
- ▶ Controlled and limited current with selected components
- ▶ The use of frequency signal splitting with duty cycles
- ▶ 5 types of signal strength in normal condition
- ▶ 10 levels with supplied cable inserted with reference disposable electrode
- ▶ Controlled output signal
- ▶ Disposable needles
- ▶ Emission mode that can be set in ten intensities with on / off times of 0.5 m sec and 10 msec with duty cycle system
- ▶ Graphic display of program control

ACCESSORIES



TECHNICAL DETAILS

Dimensions: LxAxP:
• handpiece: 30x3x190mm
• charging base (ØxA): 80x53mm
Weight: 250g
Power supply: rechargeable batteries
lithium ion 11.1 V 500 mAh
Max absorption: 30mAh (HIPLUS)

Surface treated: 1 mm²
Mode of Issue:
frequency signal splitting with graphic display
for level indication
Led: battery recharge indicator
Max peak voltage: 2100 VAC